



Vacuum creates a “perfect fit” in the new **VACOcast Fit**. It’s incredibly stable (and versatile!).

MEDICAL INDICATIONS

- Ankle fractures
- Arthrodeses of the upper and lower ankle joint
- Calcaneus fractures
- Complex ligament injuries
- Diabetic foot ulcers*
- Metatarsal fractures
- Severe sprains
- Soft tissue injuries with impaired healing process
- Splinting of foot and lower leg

*additional foam insole recommended

To order:

Contact Ken Dupuis
OPED USA
877.673.3426
ken@opedusa.com

Size Article Numbers

Small: VPCA-01-1-S
Medium: VPCA-01-1-M
Large: VPCA-01-1-L

Why vacuum?

Micro beads gently contours to the limb, as air gets sucked out of VACOcast’s patented *Vacuum Relief* insert. This creates a comfortable “perfect fit” with cast-like stability, to enable earlier weight bearing (when appropriate).

There is no inward pressure on the limb, which makes VACOcast notably comfortable—and extremely functional.

2012-02-20

How does it work? The air gets sucked out the *Vacuum Relief* insert to let micro beads contour to the leg. Because there is no inward pressure, pain is substantially minimized, enabling patients to continue being active—and heal faster.

The new VACOcast Fit is a fixed plantar-flexion boot superbly designed and engineered to help patients recover faster—yet be remarkably comfortable during their recovery.

VACOcast is also the most versatile orthopedic boot in the marketplace...

3-in-1 usability

Splint

Creates a rigid structure within seconds, to secure the limb.

Cast

Forms a sturdy cast in less than 3-minutes with no waste.

Walker

Stability and style work in unison for patients to walk in comfort.

Other features

Adjustable calf width to best fit patients with large calves, or to adapt to patients who have edema

Moveable toe plate to lengthen the plantar foot bed, to accommodate people with longer feet

Detachable sole to keep street dirt out of beds and off furniture

Velour lining is exchangeable & washable, for best hygiene practices. Additional liners are available for purchase

Security seal deters patients from opening VACOcast inappropriately, to ensure compliance (optional)

ClickDisk™ makes audible “click” noise to ensure the patient has controlled weight bearing (optional)